


pco.[®]

pco.edge series


Active-Cooled sCMOS Cameras






excelitas[®]

pco.edge series

The pco.edge series represents the high-end segment of PCO's scientific CMOS (sCMOS) camera portfolio, delivering exceptional imaging performance across a wide range of applications. Built around temperature-stabilized, high-performance sCMOS image sensors, these cameras offer extremely low readout noise, wide dynamic range, high frame rates, and high resolution across a large field of view. Across the series, selected models offer features such as low light mode, line scanning mode, double shutter mode, and lens control for specialized applications.



Technical table ¹	pco.edge 26		pco.edge 10 bi	pco.edge 10 bi LT	pco.edge 9.4 bi
Interface	CLHS FOL	USB 3.1 Gen 1	CLHS FOL		CLHS FOL
Sensor technology	sCMOS		Back illuminated sCMOS		Back illuminated sCMOS
Color type	Monochrome or color	Monochrome	Monochrome		Monochrome
Resolution [pixels]	5120 × 5120		4416 × 2368		4096 × 2300
Sensor diagonal [mm]	18.1		23.0		21.6
Pixel size [μm]	2.5 × 2.5		4.6 × 4.6		4.6 × 4.6
Max. frame rate [fps] @ full resolution	150	6	122		122
Max. pixel rate [MPixel/s]	3932	187	1467		1467
Peak quantum efficiency	65 % @ 500 nm ⁵		85 % @ 500 nm		85 % @ 500 nm
Typ. readout noise ² [e ⁻]	3.0	2.3	0.7	1.3	0.3
Dark current [e ⁻ /pixel/s] @ sensor temperature	0.4 @ +15 °C	0.09 @ -10 °C	0.4 @ +10 °C		0.4 @ +10 °C
Max. dynamic range	1300 : 1	2000 : 1	14000 : 1	14000 : 1	23300 : 1
Shutter mode ³	GS		RS		RS
Cooling method ⁴	Air & liquid		Air & liquid		Air & liquid
Features	Double shutter ⁶ , lens control, NIR-enhanced sensor ⁶	NIR-enhanced sensor ⁶	Line scanning mode, lens control	Lens control	Line scanning mode, lens control
Dimensions H x W x L [mm]	95 × 90 × 109	85 × 80 × 109	95 × 90 × 109		
Camera housing	 				

¹ Note: All values represent maximum achievable specifications in specific operational modes and may not be attainable simultaneously.

² The readout noise values are given as median (med). All values are raw data without any filtering.

³ RS = Rolling Shutter | GR = Global Reset | GS = Global Shutter



High performance through optimized mechanics

The pco.edge series features a sophisticated mechanical housing. An optimized cooling concept enables thermal stabilization of the sensor at low temperatures, resulting in a negligible dark current and improved image quality. Cooling is achieved either via external liquid cooling or through optimized internal airflow within the camera. Additionally, the cooling system's special design protects the sensor from vibrations, ensuring unimpaired imaging performance.



Technical table ¹	pco.edge 5.5		pco.edge 4.2 bi	pco.edge 4.2 bi UV
Interface	CLHS FOL	USB 3.0	USB 3.1 Gen 1	
Sensor technology	sCMOS		Back illuminated sCMOS	
Color type	Monochrome or color		Monochrome	
Resolution [pixels]	2560 × 2160		2048 × 2048	
Sensor diagonal [mm]	21.8		18.8	
Pixel size [μm]	6.5 × 6.5		6.5 × 6.5	
Max. frame rate [fps] @ full resolution	100	30	40	
Max. pixel rate [MPixel/s]	572	320	184	
Peak quantum efficiency	60 % @ 600 nm ⁵		95 % @ 580 nm	89 % @ 580 nm 48 % @ 240 nm
Typ. readout noise ² [e ⁻]	1.0		1.0	
Dark current [e ⁻ /pixel/s] @ sensor temperature	RS/GR: 0.6 GS: 0.9 @ +7 °C	RS/GR: 0.5 GS: 0.8 @ +5 °C	0.2 @ -25 °C	
Max. dynamic range	30 000 : 1		26 500 : 1	
Shutter mode ³	RS, GR, GS		RS, GR	
Cooling method ⁴	Air, optional: liquid	Air & liquid	Air & liquid	
Features	Double shutter ⁶ , lens control	-	Line scanning mode, low light mode	
Dimensions H x W x L [mm]	76 × 70 × 122.5	76 × 70 × 99.5	85 × 80 × 109	
Camera housing				

⁴Air = air forced with fan | liquid = external liquid connection

⁵ Monochrome version

⁶ Specific camera variant



Top feature – Camera Link HS



Camera Link HS (CLHS) is a standardized protocol renowned for its exceptional speed, reliability, and bandwidth. It evolved from Camera Link, the vision industry's first standard camera interface protocol. Camera Link HS is poised to provide significantly expanded capabilities in upcoming releases, positioning it as a leading choice for camera interfacing. Cameras of the pco.edge series use this mature and robust interface combined with a fiber-optic link (FOL), enabling high-speed data transmission over long distances.




Technical table ¹	pco.edge 4.2		pco.edge 4.2 LT	pco.edge 3.1
Interface	CLHS FOL	USB 3.0	USB 3.0	USB 3.0
Sensor technology	sCMOS		sCMOS	sCMOS
Color type	Monochrome		Monochrome	Monochrome or color
Resolution [pixels]	2048 × 2048		2048 × 2048	2048 × 1536
Sensor diagonal [mm]	18.8		18.8	16.6
Pixel size [μm]	6.5 × 6.5		6.5 × 6.5	6.5 × 6.5
Max. frame rate [fps] @ full resolution	100	40	40	50
Max. pixel rate [MPixel/s]	548	220	220	408
Peak quantum efficiency	80 % @ 580 nm		80 % @ 580 nm	60 % @ 600 nm ⁵
Typ. readout noise ² [e ⁻]	0.8		0.8	1.1
Dark current [e ⁻ /pixel/s] @ sensor temperature	0.6 @ +7 °C	0.3 @ 0 °C	0.8 @ +10 °C	RS/GR: 0.5 GS: 0.8 @ +5 °C
Max. dynamic range	37500 : 1		37500 : 1	27000 : 1
Shutter mode ³	RS	RS, GR	RS, GR	RS, GR, GS
Cooling method ⁴	Air, optional: liquid	Air & liquid	Air	Air
Features	Lens control	-	-	-
Dimensions H x W x L [mm]	76 × 70 × 122.5	76 × 70 × 99.5	76 × 70 × 99.5	76 × 70 × 99.5
Camera housing				



Top feature – Line scanning mode


The PCO line scanning mode is a special readout mode dedicated to lightsheet fluorescence microscopy (LSFM), enabling optimized synchronization of the camera and the lightsheet microscope system. This feature is built on the rolling shutter mode. Compared to the standard rolling shutter, in line scanning mode parameters for the number of exposure lines and line time are adjustable. The number of exposure lines corresponds to the slit width, while the line time defines the slit speed.

The camera supports various trigger options for external synchronization.




Exemplary illustration of the readout in line scanning mode with five exposure lines.

Quantum efficiency




Top feature – Low light mode

In standard mode, two images with exactly the same exposure time but different gains are recorded. The low gain image is optimized for high fullwell capacity and the high gain image is optimized for low readout noise. Both images are merged into one high dynamic range image.



The low light mode benefits from two times correlated multi sampling of high gain images. This reduces the temporal noise by a factor of the square root of 2, which is ideal for applications demanding low noise and high sensitivity.



About Excelitas®

Excelitas is a leading provider of advanced, life-enriching technologies that make a difference, serving global market leaders in the life sciences, advanced industrial, next-generation semiconductor and avionics end markets. Headquartered in Pittsburgh, PA, USA, Excelitas is an essential partner in the design, development and manufacture of advanced technologies, offering leading-edge innovation in sensing, detection, imaging, optics and specialty illumination for customers worldwide.

Excelitas is at the forefront of addressing many of the relevant megatrends impacting the world today, including precision medicine, industrial automation, artificial intelligence, and connected devices (IoT).



+49 9441 2005 0
Europe

+1 866 662 6653
North America

+86 400 080 2255
Asia-Pacific

Donaupark 11
93309 Kelheim
Germany

excelitas.com
pco@excelitas.com

For a complete listing of our global offices, visit www.excelitas.com/locations

© 2026 Excelitas Technologies Corp. All rights reserved. The Excelitas logo, Excelitas® and PCO® are registered trademarks of the Excelitas group of companies. All other products and services are either trademarks or registered trademarks of their respective owners. Excelitas reserves the right to change this document at any time without notice and disclaims liability for editorial, pictorial or typographical errors.

PO_pco.edge_v2.0.0