## SPCM-AQRH-XX-TR

# **Timing Resolution Optimized Single Photon Counting Module**



Excelitas Technologies' SPCM-AQRH-XX-TR is a Single Photon Counting Module of the most recent product generation, specifically selected and performanceoptimized for timing resolution.

The SPCM-AQRH-XX-TR uses a specially selected SLiK silicon avalanche photodiode giving module timing resolution better than 250 ps while maintaining peak photon detection efficiency (PDE) of more than 75% at 650 nm over a 180  $\mu$ m diameter active area. While some performance parameters, in particular afterpulse probability, are traded off against the optimized timing performance, other performance parameters of the standard SPCM-AQRH, such as outstanding uniformity, overload protection, temperature stability and linearity, are still maintained by this new timing-optimized module.

This family of –TR modules is designed to support applications in time correlated single photon counting (TCSPC), fluorescence lifetime measurements and fluorescence lifetime imaging microscopy (FLIM).

Excelitas' series of photon counting modules are designed and built to be fully compliant with the European Union Directive 2015/863/EU – Restriction of the use of certain Hazardous Substances in Electrical and Electronic equipment (RoHS).

#### **Key Features**

- Timing resolution <250ps</li>
- Peak photon detection efficiency (PDE) @ 650 nm: 75% typical
- Active area: 180 μm
- Gated output
- Single +5 V supply
- RoHS-compliant
- Linearity over high count rate

#### **Applications**

- Time correlated single photon counting
- Fluorescence lifetime imaging microscopy
- Ultra-sensitive fluorescence lifetime measurements
- Quantum Cryptography
- Photon correlation spectroscopy
- Optical range finding
- Particle sizing
- Adaptive Optics

## **SPCM-AQRH-XX-TR Series**

# **Timing Resolution Optimized Single Photon Counting Module**

Table 1. Specifications of SPCM-AQRH-XX-TR, @ 22 °C, all models; unless otherwise indicated (1)

Parameter		<b>Min</b> 170	Тур	Max	<b>Unit</b> μm	
Active area (dia	meter) at minimum PDE		180			
Photon detection	on efficiency (PDE) (2) at:					
	650 nm		75		%	
	830 nm		50		%	
Dark Count	SPCM-AQRH-W0			1500		
	SPCM-AQRH-W1			1000		
	SPCM-AQRH-W2			500	Counts /	
	SPCM-AQRH-W3			250	second	
	SPCM-AQRH-W4			100		
Single photon t	iming resolution (at 825 nm) <sup>(2,3)</sup>					
	t factory for optimized timing below		225	250	ps	
200 ps and at other wavelengths					·	
Output pulse w	idth <sup>(9)</sup>					
	SPCM-AQRH-1X, SPCM-AQRH-4X		10		ns	
	SPCM-AQRH-2X, SPCM-AQRH-5X		18		ns	
	SPCM-AQRH-3X, SPCM-AQRH-6X		28		ns	
See table 3.						
Dead time (cou	nt rate below 5M/c)					
	SPCM-AQRH-1X, SPCM-AQRH-4X		22		ns	
	SPCM-AQRH-2X, SPCM-AQRH-5X		28		ns	
	SPCM-AQRH-3X, SPCM-AQRH-6X		35		ns	
See table 3.						
Output pulse ar	nplitude:					
SPCM-AQRH-1>	X, SPCM-AQRH-2X, SPCM-AQRH-3X					
	TTL HIGH	1.5	2.2		V	
	TTL LOW	-0.1		0.8	V	
See table 3.						
SPCM-AQRH-4X	K, SPCM-AQRH-5X, SPCM-AQRH-6X					
	TTL HIGH	3.0	4.4		V	
	TTL LOW	-0.1	7.4	0.8	V	
See table 3.		-0.1		0.8	V	
Linearity correc			1			
	1 Mc/s		1.02			
	5 Mc/s		1.16			
	10 Mc/s		1.40			
Afterpulsing pro	obability		1.0	3.0	%	

<sup>(1)</sup> For other performance characteristics, refer to Operating Instructions, product notes and specifications listed on the standard SPCM-AQRH data sheet.

<sup>(2)</sup> See Figure 1 for timing resolution vs. photon detection efficiency curve.

<sup>(3)</sup> Timing resolution is measured using a 10 µm diameter light spot, at 825 nm, case temperature at 22°C. For timing resolution requirements measured with a larger spot size, or at different wavelength, please contact Excelitas.

## **SPCM-AQRH-XX-TR Series**

# **Timing Resolution Optimized Single Photon Counting Module**

## **Table 2. Absolute Maximum Ratings**

Supply voltage (1)	5.5 V
Maximum count rate	Maximum count rate can be sustained if case temperature is maintained within specified limits.
Peak light intensity	10⁴ photons per pulse and pulse width < 1 ns
Case temperature (3)	-20°C/+70°C storage, +5°C /+70°C operating

**Table 3. SPCM Ordering Guide** 

Order Part# W - Output Pulse Options			X - Dark Count Rates					
WX-YY	Output Pulse	Dead Time	Output Pulse	-W0	-W1	-W2	-W3	-W4
	Width (ns)	(ns)	Height (V)					
SPCM-AQRH-1X-TR	10	22	2.2	≤1500 cps	≤1000 cps	≤500 cps	≤250 cps	≤100 cps
SPCM-AQRH-2X-TR	18	28	2.2	≤1500 cps	≤1000 cps	≤500 cps	≤250 cps	≤100 cps
SPCM-AQRH-3X-TR	28	35	2.2	≤1500 cps	≤1000 cps	≤500 cps	≤250 cps	≤100 cps
SPCM-AQRH-4X-TR	10	22	4.4	≤1500 cps	≤1000 cps	≤500 cps	≤250 cps	≤100 cps
SPCM-AQRH-5X-TR	18	28	4.4	≤1500 cps	≤1000 cps	≤500 cps	≤250 cps	≤100 cps
SPCM-AQRH-6X-TR	28	35	4.4	≤1500 cps	≤1000 cps	≤500 cps	≤250 cps	≤100 cps

Example: SPCM-AQRH-43-TR = 10ns output pulse width, 22ns dead time, 4.4V output pulse height, <250 cps dark count rate




Figure 1: Typical timing resolution & PDE relative to position of detector chip @ 650 nm

## **Timing Resolution Optimized Single Photon Counting Module**

## Warranty

A standard 12-month warranty following shipment applies. Any warranty is null and void if the module case has been opened. Warranty is null and void if the module input exceeds 5.5V or the polarity of the +5V supply is reversed.

### **Individual Module Test Data**

Each module is supplied with test data indicating the module's actual dark count, dead time, pulse width, photon detection efficiency @ 630 nm, timing resolution, and linearity correction factor.

## **About Excelitas Technologies**

Excelitas Technologies is a global technology leader focused on delivering innovative, customized solutions to meet the detection, lighting, and other high-performance technology needs of OEM customers.

From analytical instrumentation to clinical diagnostics, medical, industrial, safety and security, and aerospace and defense applications, Excelitas Technologies is committed to enabling our customers' success in their end-markets. Excelitas Technologies has approximately 3,000 employees in North America, Europe and Asia, serving customers across the world.

**Excelitas Technologies** 

22001 Dumberry Road Vaudreuil-Dorion, Quebec Canada J7V 8P7 Telephone: (+1) 450.424.3300 Toll-free: (+1) 800.775.6786 Fax: (+1) 450.424.3345 detection.na@excelitas.com Excelitas Technologies GmbH & Co. KG Wenzel-Jaksch-Str. 31

D-65199 Wiesbaden Germany Telephone: (+49) 611 492 430 Fax: (+49) 611 492 165

Telephone: (+65) 6775 2022 (Main number)

Fax: (+49) 611 492 165 (Customer Service) detection.europe@excelitas.com Fax: (+65) 6778-1752

Fax: (+65) 6778-1752 detection.asia@excelitas.com

Telephone: (+65) 6770 4366

**Excelitas Technologies** 

8 Tractor Road

Singapore 627969



For a complete listing of our global offices, visit <a href="www.excelitas.com/Locations">www.excelitas.com/Locations</a>

© 2020 Excelitas Technologies Corp. All rights reserved. The Excelitas logo and design are registered trademarks of Excelitas Technologies Corp. All other trademarks not owned by Excelitas Technologies or its subsidiaries that are depicted herein are the property of their respective owners. Excelitas reserves the right to change this document at any time without notice and disclaims liability for editorial, pictorial or typographical errors.