

X-Cite® illuminators have been qualified for the following microscopes:

Microscope Model	Light Guide Adaptor (for X-Cite 120Q series, 200DC, exacte, XLED1, 110LED, TURBO)	Direct Coupling Flange (for X-Cite 120LED/Boost /mini)
CARL ZEISS		
AxioExaminer	810-00022X	810-10022X
AxioImager	810-00022X	810-10022X
AxioObserver	810-00022X	810-10022X
AxioScope A1	810-00022X	810-10022X
Axiophot, Axiophot2	810-00022X	810-10022X
Axioplan 1, 2	810-00022X	810-10022X
Axioskop 1, 20, 50	810-00022X	810-10022X
Axioskop 2, 2FS, 40	810-00022X / 810-00035X†	810-10022X
AxioTech	810-00022X	810-10022X
Axiovert A1	810-00022X	810-10022X
Axiovert 10, 35, 100, 135, S100TV	810-00035X	TBD
Axiovert 200	810-00022X / 810-00035X†	810-10022X
Axiovert 400	810-00035X	TBD
Axiovert 25, 40	810-00037	810-10022X
AxioZoom	No adaptor req'd*	N/A
Discovery.v8 / .v12 / .v20	No adaptor req'd*	N/A
LSM	No adaptor req'd*	N/A
Stemi SV11	No adaptor req'd*	N/A
Stereo Lumar	810-00022X	810-10022X
LEICA		
DMi8 – Manual	810-00033X	810-10033X
DMi8 – Automated w/FIM	No adaptor req'd*	N/A
DM2000, DM6000, DMI4000	810-00033X	810-10033X
DMIL	810-00033X	810-10033X
DMIRB, DMIRM, DMIRE2	810-00033X	810-10033X
DMLB	810-00033X	810-10033X
DMLFS	810-00033X	810-10033X
DMR, DMRXA, DMRXA2	810-00033X	810-10033X
MZFLIII, MZ16F, MZ16FA	810-00038X	810-10038X
M165FC	810-00038X	810-10038X
M10, MZ9.5, MZ10F, M205FA	810-00038X	810-10038X
MOTIC		
AE31	810-00055	N/A
BA400	810-00055	N/A

NIKON		
50i, 80i, 90i	810-00030X	810-10030X
AZ100	No adaptor req'd*	N/A
Diaphot, Labophot, TMF	810-00030X*	TBD
E400, E600, E800	810-00030X	810-10030X
FN-1	810-00030X	810-10030X
SMZ1500, SMZ800	810-00030X	810-10030X
SMZ18, SMZ25	No adaptor req'd*	N/A
TE2000, all other "TE"	810-00030X	810-10030X
Ti, Ni, Ci	810-00030X	810-10030X
TS100	810-00030X*	810-10093X
For "white light TIRF"	810-00045	TBD
most other Nikon microscopes	810-00030X	810-10030X
OLYMPUS		
AX70	810-00032X	TBD
BX series: BX50, BX60	810-00032X	810-10032X
BX2 series: BX41, BX51, BX61	810-00032X	810-10032X
BX3 series: BX43, BX46, BX53	810-00032X	810-10032X
CKX41	810-00032X	810-10032X
CK40	810-00032X	TBD
CX41	810-00077X	TBD
IX series: IX50, IX70	810-00031X	810-10031X
IX2 series: IX51, IX71, IX81	810-00031X	810-10031X
IX3 series (w/ RFA)	810-00046X	810-10031X
IX3 series (w/ RFA-L)	810-00046X	810-10032X
MVX10	810-00046X	810-10031X
SZX7, 10, 12, 16	810-00032X	810-10032X

Notes:

* Additional fittings and/or specific configuration from microscope manufacturer may be required.

† With external devices (e.g. shutters, filter wheels, etc) installed between the adaptor and microscope, use 810-00035X.