

X-Cite® illuminators have been qualified for the following microscopes:

Microscope Model	Light Guide Adaptor For: X-Cite 120Q series, XLED1, 110LED, XYLIS II, NOVEM Also: 200DC, exacte, FIRE, TURBO, XYLIS	Direct Coupling Flange For: X-Cite mini+ Also: 120LED, 120LEDmini	Direct Coupling Flange For: X-Cite 120LEDBoost
CARL ZEISS			
AxioExaminer	810-00022X	810-10022X	810-10022X
Axiomager	810-00022X	810-10022X	810-10022X
AxioObserver	810-00022X	810-10022X	810-10022X
AxioScope A1	810-00022X	810-10022X	810-10022X
AxioScope 5	810-00022X	810-10022X	810-10022X
AxioScan	810-00022X	810-10022X	Not recommended
AxioVert A1	810-00022X	810-10022X	810-10022X
AxioZoom	No adaptor req'd*	N/A	N/A
Discovery.v8 / .v12 / .v20	No adaptor req'd*	N/A	N/A
LSM	No adaptor req'd*	N/A	N/A
M2Bio (w/Kramer attachment)	810-00022X	810-10022X	TBD
<i>Axiophot, Axiophot2</i>	<i>810-00022X</i>	<i>810-10022X</i>	<i>810-10022X</i>
<i>Axioplan 1, 2</i>	<i>810-00022X</i>	<i>810-10022X</i>	<i>810-10022X</i>
<i>Axioskop 1, 20, 50</i>	<i>810-00022X</i>	<i>810-10022X</i>	<i>810-10022X</i>
<i>Axioskop 2, 2FS, 40</i>	<i>810-00022X / 810-00035X†</i>	<i>810-10022X</i>	<i>810-10022X</i>
<i>AxioTech</i>	<i>810-00022X</i>	<i>810-10022X</i>	<i>TBD</i>
<i>Axiovert 10, 35, 100, 135, S100TV</i>	<i>810-00035X</i>	<i>TBD</i>	<i>TBD</i>
<i>Axiovert 200</i>	<i>810-00022X / 810-00035X†</i>	<i>810-10022X</i>	<i>TBD</i>
<i>Axiovert 25, 40</i>	<i>810-00037</i>	<i>810-10022X</i>	<i>TBD</i>
<i>Stemi SV11</i>	<i>No adaptor req'd*</i>	<i>N/A</i>	<i>N/A</i>
<i>Stereo Lumar</i>	<i>810-00022X</i>	<i>810-10022X</i>	<i>TBD</i>

Notes:

* Additional fittings and/or specific configuration from microscope manufacturer may be required.

† With external devices (e.g. shutters, filter wheels, etc) installed between the adaptor and microscope, use 810-00035X.

X-Cite and microscope models in italics have been discontinued.

X-Cite® illuminators have been qualified for the following microscopes:

Microscope Model	Light Guide Adaptor For: X-Cite 120Q series, XLED1, 110LED, XYLIS II, NOVEM Also: 200DC, exacte, FIRE, TURBO, XYLIS	Direct Coupling Flange For: X-Cite mini+ Also: 120LED, 120LEDmini	Direct Coupling Flange For: X-Cite 120LEDBoost
LEICA			
DM4B, DM6B	810-00033X	Not recommended	810-10033X + 8mm spacer
DMi8 – Manual	810-00033X	810-10033X	TBD
DMi8 – Automated w/FIM	No adaptor req'd*	N/A	N/A
DMIL	810-00033X	TBD	TBD
DM2000/2500/3000	810-00033X	TBD	TBD
M165FC, M205FA	810-00038X	810-10038X‡	TBD
MZ10F	810-00038X	TBD	TBD

<i>DM4000-6000</i>	<i>810-00033X</i>	<i>TBD</i>	<i>TBD</i>
<i>DMI4000</i>	<i>810-00033X</i>	<i>TBD</i>	<i>TBD</i>
<i>DMI6000B</i>	<i>810-00033X</i>	<i>810-10033X</i>	<i>TBD</i>
<i>DMIRB, DMIRM, DMIRE2</i>	<i>810-00033X</i>	<i>Not recommended</i>	<i>TBD</i>
<i>DMLB</i>	<i>810-00033X</i>	<i>TBD</i>	<i>TBD</i>
<i>DMLFS</i>	<i>810-00033X</i>	<i>TBD</i>	<i>TBD</i>
<i>DMR, DMRXA, DMRXA2</i>	<i>810-00033X</i>	<i>TBD</i>	<i>TBD</i>
<i>MZFLIII</i>	<i>810-00038X</i>	<i>810-10038X‡</i>	<i>TBD</i>
<i>MZ16F, MZ16FA</i>	<i>810-00038X</i>	<i>TBD</i>	<i>TBD</i>
<i>M10, MZ9.5</i>	<i>810-00038X</i>	<i>TBD</i>	<i>TBD</i>
<i>Leitz - all models</i>	<i>Incompatible</i>	<i>Incompatible</i>	<i>Incompatible</i>

Notes:

* Additional fittings and/or specific configuration from microscope manufacturer may be required.

‡ May not be suitable for some applications, demo before purchase recommended.

X-Cite and microscope models in italics have been discontinued.

X-Cite® illuminators have been qualified for the following microscopes:

Microscope Model	Light Guide Adaptor For: X-Cite 120Q series, XLED1, 110LED, XYLIS II, NOVEM Also: 200DC, exacte, FIRE, TURBO, XYLIS	Direct Coupling Flange For: X-Cite mini+ Also: 120LED, 120LEDmini	Direct Coupling Flange For: X-Cite 120LEDBoost
NIKON			
AZ100	Requires Nikon Part: MNE84010	N/A	N/A
Ci (w/HMX connector)	810-00030X	810-10030X	Not recommended
Ci (w/MBE94350 EPI-FL attachment)	Requires Nikon Part: MXA22205	TBD	
FN1 (w/HMX connector)	810-00030X	810-10030X	TBD
LV100	N/A	810-10097X	TBD
Ni (w/HMX connector)	810-00030X	810-10030X	Not recommended
Ti (w/HMX connector)	810-00030X	810-10030X	Not recommended
Ti2 (w/Ti-LAPP)	No adaptor req'd*	N/A	N/A
Ti2 (w/MEE54740 EPI-FL attachment)	Requires Nikon Part: MXA22205	TBD	
SMZ18, SMZ25	No adaptor req'd*	N/A	N/A
SMZ1270i/1270/800N	No adaptor req'd*	N/A	N/A
For "white light TIRF"	810-00045	N/A	N/A
<i>50i, 80i, 90i</i>	<i>810-00030X</i>	<i>810-10030X</i>	<i>TBD</i>
<i>Diaphot, Labophot, TMF</i>	<i>810-00030X*</i>	<i>810-10030X*</i>	<i>TBD</i>
<i>E400, E600, E800</i>	<i>810-00030X</i>	<i>810-10030X</i>	<i>Not recommended</i>
<i>SMZ800, SMZ1500</i>	<i>810-00030X</i>	<i>810-10030X</i>	<i>TBD</i>
<i>TE2000, all other "TE"</i>	<i>810-00030X</i>	<i>810-10030X</i>	<i>TBD</i>
<i>TS100</i>	<i>810-00030X*</i>	<i>810-10093X</i>	<i>TBD</i>

Notes:

* Additional fittings and/or specific configuration from microscope manufacturer may be required.

‡ May not be suitable for some applications, demo before purchase recommended.

X-Cite and microscope models in italics have been discontinued.

X-Cite® illuminators have been qualified for the following microscopes:

Microscope Model	Light Guide Adaptor For: X-Cite 120Q series, XLED1, 110LED, XYLIS II, NOVEM Also: 200DC, exacte, FIRE, TURBO, XYLIS	Direct Coupling Flange For: X-Cite mini+ Also: 120LED, 120LEDmini	Direct Coupling Flange For: X-Cite 120LEDBoost
OLYMPUS / EVIDENT			
BX2 series: BX51WI, BX61WI	810-00032X	810-10032X	810-10032X
BX3 series: BX43, BX46, BX53, BX63	810-00032X	810-10032X	810-10032X
CKX53	810-00032X	810-10032X	TBD
IX3 series (w/ RFA)	810-00046X	810-10031X	Not recommended
IX3 series (w/ RFA-L)	810-00046X	810-10032X	810-10032X-01
MVX10	810-00046X	810-10031X	TBD
SZX7, 10, 16	810-00032X	810-10032X	TBD
VS200	810-00046X	N/A	N/A

<i>AX70</i>	<i>810-00032X</i>	<i>TBD</i>	<i>TBD</i>
<i>BX series: BX50, BX60</i>	<i>810-00032X</i>	<i>810-10032X</i>	<i>810-10032X</i>
<i>BX2 series: BX41, BX51, BX61</i>	<i>810-00032X</i>	<i>810-10032X</i>	<i>810-10032X</i>
<i>CKX41</i>	<i>810-00032X</i>	<i>810-10032X</i>	<i>TBD</i>
<i>CK40</i>	<i>810-00032X</i>	<i>TBD</i>	<i>TBD</i>
<i>CX41</i>	<i>810-00077X</i>	<i>TBD</i>	<i>TBD</i>
<i>IX series: IX50, IX70</i>	<i>810-00031X</i>	<i>810-10031X</i>	<i>TBD</i>
<i>IX2 series: IX51, IX71, IX81</i>	<i>810-00031X</i>	<i>810-10031X</i>	<i>TBD</i>
<i>SZX12</i>	<i>810-00032X</i>	<i>810-10032X</i>	<i>TBD</i>
<i>VS120</i>	<i>810-00032X</i>	<i>N/A</i>	<i>N/A</i>

Notes:

X-Cite and microscope models in italics have been discontinued.

X-Cite® illuminators have been qualified for the following microscopes:

Microscope Model	Light Guide Adaptor For: X-Cite 120Q series, XLED1, 110LED, XYLIS II, NOVEM <i>Also: 200DC, exacte, FIRE, TURBO, XYLIS</i>	Direct Coupling Flange For: X-Cite mini+ <i>Also: 120LED, 120LEDmini</i>	Direct Coupling Flange <i>For: X-Cite 120LEDBoost</i>
ACCU-SCOPE			
EXC-400	810-00032X	810-10096X	TBD
EXC-500	810-00032X	810-10032X	TBD
EXI-310	810-00032X	810-10032X	TBD
EXI-600	810-00022X	810-10022X	TBD

Notes:

X-Cite and microscope models in italics have been discontinued.