

Filter Recommendations

Fluorophore	Filter Type	Semrock Filter	Chroma Filter
DAPI	Single	LED-DAPI-A	49028
CFP	Single	LED-CFP-A	49001
CFP/YFP FRET	FRET multiple set	FRET CFP/YFP-C	89002
FITC or GFP	Single	LED-FITC-A	49002
YFP	Single	LED-YFP-A	49003
TRITC	Single	LED-TRITC-A	49004
mCherry	Single	LED-mCherry-A	49008
Cy5	Single	LED-Cy5-A	49006