

## X-Cite® 120LEDmini Filter Recommendations

Fluorophore	Filter Type	Semrock Filter	Chroma Filter
DAPI	Single	<a href="#">LED-DAPI-A</a>	49028
CFP	Single	<a href="#">LED-CFP-A</a>	49001
CFP/YFP FRET	FRET multiple set	<a href="#">FRET CFP/YFP-C</a>	89002
FITC or GFP	Single	<a href="#">LED-FITC-A</a>	49002
YFP	Single	<a href="#">LED-YFP-A</a>	49003
TRITC	Single	<a href="#">LED-TRITC-A</a>	49004
mCherry	Single	<a href="#">LED-mCherry-A</a>	49008
Cy5	Single	<a href="#">LED-Cy5-A</a>	49006