

To simplify system selection, the XLED1 offers three bundled configurations as standard options that are designed to address primary fluorescence applications.

LED Module	Excitation Range of Use	Bundled System Configurations		
		Live Cell	Fixed Cell #1	Fixed Cell #2
UVA	360-380 nm		•	
UVX	375-400 nm			•
UVV	390-415 nm			
BLX	430-465 nm	•		
BDX	450-495 nm	•	•	•
BGX	505-545 nm	•		
GYX	540-600 nm	•	•	•
RLX	615-655 nm		•	•
RDX	660-675 nm			
IRX	710-750 nm			