

Software Supporting X-Cite Products

Manufacturer	Software	X-Cite 120PC Q	X-Cite 200DC	X-Cite exacte	X-Cite XLED1	X-Cite 110LED	X-Cite 120LED Series / mini+	X-Cite XYLIS	X-Cite TURBO	X-Cite NOVEM
Carl Zeiss	Zen - Black	•								
	Zen - Blue			•	○	•	•	•	○	○
Explora Nova	Fluo Up Mercator Morpho Strider	•		•		C	C	C		
Hamamatsu	HCIImage	•		•	•	•	•	•	•	•
Intelligent Imaging Innovations (3i)	Slidebook	•	○	•	•	*	*	*	+	•
Leica	LAS-X		○	○	•	○	○	#	+	
	Cytovision	•		*		*	*	C		
Media Cybernetics	Scope-Pro / Image-Pro	•		•	•	C	C	C	+	
MicroBrightField (MBF)	NeuroLucida Stereo Investigator MicroLucida	•		•	•	L	•	•	•	•
MicroManager	MicroManager	•		•	•	*	*	*	•	•
Molecular Devices	MetaMorph	•	○	•	•	*	*	C	•	•
Nikon	NIS - Elements	•	○	•	•	*	*	•	•	T*
Olympus	cellSens	•	○	•	•	•	•	•	•	T
	VS-ASW for VS120			•		*	*			
	VS-ASW for VS200							•	•	•

Jan 2022
MG

- Driver Available
- TTL Shutter only
- * Compatible with X-Cite 120PC/exacte driver (note: if available, use exacte driver for higher level of control)
- C Compatibility with X-Cite 120PC/exacte driver expected
- + Compatible with XLED1 driver for up to 4 channels, Grouping Mode available
- # Compatible with XLED1 driver - Channel 1
- L Compatible with 120LED driver
- T Compatible with TURBO driver for up to 6 channels, Grouping Mode available
- T* Compatible with TURBO driver for channels 1-6 only, Grouping Mode not available