

## X-Cite® XYLIS Control Options

### Control Options

|                            |   |
|----------------------------|---|
| Manual Control – speedDIAL | ON/OFF<br>Speed sensitive intensity control<br>1 favorite intensity setting<br>UV OFF<br>Left and right-handed control option |
| USB                        | ON/OFF<br>Intensity control<br>Compatible with X-Cite 120PC, X-Cite exacte and X-Cite 120LED drivers                          |
| TTL                        | ON/OFF  |
| Foot Pedal (optional)      | ON/OFF<br>All X-Cite XYLIS models have input port, foot pedal sold separately   |