


data sheet

pco. lens adapter EF-mount

This **lens control adapter** enables the user to connect electronic EF- and EF-S / Canon mount lenses to a variety of PCO cameras. The adapter allows remote control focus and aperture of compatible lenses.

Use our **MachVis** lens selector software to learn more about our cameras and their optical interface.



compatible PCO cameras

pco.edge


pco.edge 26 CLHS, pco.edge 10 bi CLHS, pco.edge 10 bi LT CLHS, pco.edge 5.5 CLHS, pco.edge 5.5 DS CLHS, pco.edge 4.2 CLHS

pco.dicam

pco.dicam C1, pco.dicam C1 LX, pco.dicam C1 LT, pco.dicam C1 UHS, pco.dicam C4, pco.dicam C4 LX, pco.dicam C4 LT, pco.dicam C4 UHS, pco.dicam C8, pco.dicam C8 LX, pco.dicam C8 LT, pco.dicam C8 UHS

pco.dimax


pco.dimax cs4, pco.dimax cs3, pco.dimax cs1



pin assignment

pin 1	Power input	12 V / 15 W
pin 2	Signal output	TXD, communication for lens control
pin 3	Signal input	RXD, communication for lens control
pin 4	Ground	

dimensions



ordering information

pco.lens adapter EF-mount

30108025020

EF-mount adapter with lens control function; compatible with pco.edge, pco.dimax and pco.dicam series