

MEDIA ADVISORY

November 20, 2019

**Excelitas to Display OmniCure LED UV Curing Systems for
Medical Device Applications at BIOMEDevice**

WHO: [Excelitas Technologies®](#), a global technology leader delivering innovative, customized photonic solutions, will showcase its innovative OmniCure® LED UV Curing Systems for medical device manufacturers at [BIOMEDevice 2019](#).

WHAT: Featured products at Excelitas booth #1043 include:

- [OmniCure S2000 Spot UV Curing Lamp System](#): The OmniCure S2000 Spot UV Curing System is ideal for precision adhesive curing in a variety of industries, particularly medical device bonding applications. The OmniCure S2000 provides fast and strong cures of UV adhesives, which result in superior product quality, rapid production and lower manufacturing costs.
- [OmniCure LX500 LED Spot UV Curing System](#): Featuring up to four LED heads with up to 16W/cm² peak irradiance and proprietary feedback technology, the LX500 provides unmatched optical stability via Intelli-Lamp® LED technology for consistent and repeatable assembly of medical devices with lower costs. OmniCure LX500 is ideal for use in medical devices such as catheters, cannulas, endoscopes and syringes.
- [OmniCure AC8-HD \(High Dose\) LED UV Curing System](#): Offering higher optical power and enhanced performance in the same air-cooled mechanical enclosure as its AC8 and AC9 predecessors, the AC8-HD System delivers over 50% more dose with unprecedented optical uniformity. Excelitas' patented control technology for addressing individual LED UV modules enables an outstandingly tight irradiation across the UV emitting length, as well as between heads. These light sources are ideal for large area curing of a variety of medical applications.

Company experts will be on site or available by phone in advance of the conference to provide updates on its solutions.

WHEN: December 4 – 5, 2019

WHERE: [San Jose Convention Center](#), San Jose, CA
Booth #1043

###

About Excelitas Technologies

Excelitas Technologies® Corp. is a photonics technology leader focused on delivering innovative, high-performance, market-driven solutions to meet the lighting, optronics, detection and optical technology needs of our OEM customers. Serving a vast array of applications across biomedical, scientific, safety, security, consumer products, semiconductor, industrial manufacturing, defense and aerospace sectors, Excelitas stands committed to enabling our customers' success in their end-markets. Our photonics team consists of 6,700 professionals



TECHNOLOGIES®

working across North America, Europe and Asia, to serve customers worldwide. Connect with Excelitas on [Facebook](#), [LinkedIn](#) and [Twitter](#).

Visit <http://www.excelitas.com> for more information.

Excelitas®, OmniCure® and Intelli-Lamp® are registered trademarks of Excelitas Technologies Corp. All other products and services are either trademarks or registered trademarks of their respective owners.

Contacts:

Scott Orr

Senior Director of Global Marketing - Commercial

scott.orr@excelitas.com

(+1) 781-996-5925

Cheryl Reynhout or Jill Anderson

On Behalf of Excelitas Technologies Corp.

SVM Public Relations

excelitas@svmmarcom.com

(+1) 401-490-9700